

Hygienic Hi-Flex Wall Bracket System, 16.5", Lime



This hygienic, modular wall bracket provides a flexible solution for storing your cleaning tools and utensils in any area of your site. The rail easily attaches to the four wall mounts by sliding it on from the right or left side. The grip band modules and hook modules can then be clicked on the rail and locked in place using one of the many small notches in the rail. These ensure that the modules cannot slide from side to side. From one to 10 hanging modules can be placed on the rail (up to five grip band modules or up to 10 hook modules). The grip band module can hold products with a diameter of 28-33 mm. The hook module can hold products weighing up to 3 kg. The system comes with three grip band modules and two hooks. Modules and rails are easy to disassemble for cleaning or replacement.

Technical Data

Item Number	101177
Holds products between	25 - 34 mm
Material	Polypropylene TPE Rubber Polyamide
Complies with (EC) 1935/2004 on food contact materials ¹	Yes
Produced according to EU Regulation 2023/2006/EC of Good Manufacturing Practice	Yes
FDA-compliant raw material (CFR 21)	Yes
Complies with UK 2019 No. 704 on food contact materials	Yes
Meets the REACH Regulation (EC) No. 1907/2006	Yes
Use of phthalates and bisphenol A	No
Is Halal and Kosher compliant	Yes
Design Registration No.	EU 008195606-0001-12
Box Quantity	5 Pcs.
Quantity per Pallet (80 x 120 x 200 cm)	360 Pcs.
Quantity Per Layer (Pallet)	30 Pcs.
Length	16.5 "
Width	3.4 "
Height	1.9 "
Net Weight	0.6349 lbs
Weight bag (Recycling Symbol "4" Low Density Polyethylene (LDPE)	0.05 lbs
Tare Cardboard	0.1 lbs
Tare Total	0.15 lbs
Gross Weight	0.7866 lbs
Cubic Feet	0.062 ft3
Recommended sterilisation temperature (Autoclave)	249.8 °F
Max. cleaning temperature (Dishwasher)	199.4 °F
Max usage temperature (food contact)	104 °F
Max usage temperature (non food contact)	176 °F
Min. usage temperature	32 °F
Min. pH-value in usage concentration	2 pH
Max. pH-value in Usage Concentration	10.5 pH
Recycling Symbol "5", Polypropylene (PP)	No
GTIN-12 Number	5705022028654
GTIN-12 Number (Box quantity)	15705028028662
Customs Tariff No.	39241000

New equipment should be cleaned, disinfected, sterilized, and any labels removed, as appropriate to its intended use, e.g. high risk vs. low risk food production areas, general hospital areas vs. intensive care units, before use.

1. See Declaration of Compliance for further details on food contact
3. Do not store the product below 32 °Fahrenheit.