

KOLD KEEPER®**FEATURES AND BENEFITS:**

- Convenient source to provide cold transport.
- Outer body of non-corrosive aluminum permanently sealed to safely hold contents.
- Standard size allows Keepers to fit any cabinet designed to hold cold foods.
- Keepers utilize non-toxic ingredients that are self-contained within outer body.
- Kold Keepers provide hours of use when pre-frozen for 8-12 hours.
- Recommended for use in cabinet series 100, 101, 102, 103, 108, 109, 150, 300, 302 and 309.




K-1218

**ACCESSORIES and OPTIONS**
(Available at extra cost):

- ☐ Keeper Carrier (Carries two Keepers)



Accessory Carrier #0645-003



CRESCOR MODEL NO.	PAN SIZE INCHES	DIM "A"		DIM "B"		DIM "C"		ACT. WGT.	
		INCHES	MM	INCHES	MM	INCHES	MM	LBS.	KG.
K-1218	12 X 20	12-7/8	330	17-7/8	455	2	55	14	6

Two will fit in the same space as one 18" x 26" pan.

- Outer body: .063 aluminum top and bottom; hemmed and crimped.
- Interior: Filled with reusable cold gel packaged in non-leaking 4 mil bag.
- Handle: Plastic covered 3/16 dia. aluminum, customer installed.
- Keepers provide 3 to 4 hours storage.
- We recommend...
 - ...4 to 6 Keepers for insulated cabinets 69-3/4" high.
 - ...6 to 8 Keepers for non-insulated cabinets 57-3/4" and 69-3/4" high.
 - ...2 to 4 Keepers for half-size and hand-lift cabinets.

SHORT FORM SPECIFICATIONS

Cres Cor Kold Keeper Model K-1218. Body formed of .063 aluminum; hemmed and crimped. Internal gel refrigerant sealed in 4 mil bag, non-toxic. Reusable, frozen cold source. NSF listed.



5925 Heisley Road • Mentor, OH 44060-1833
Phone: 877/CRESCOR • Fax: 440/350-7267
www.crescor.com

© Crescent Metal Products, Inc. 2008 All rights reserved.

In line with its policy to continually improve its products, CRES COR reserves the right to change materials and specifications without notice.

Litho in U.S.A.