Northwoods Chili

Product Specifications

Code: 63785 **Pack:** 4 x 4 pound bags (16 lbs.)

Date Effective: 6/13/18 **UPC:** 07072363785-2 **GTIN:** 0007072363785-2

Description:

A hearty chili made with premium fresh ground beef, kidney beans, onions, tomatoes, green peppers and a hint of brown sugar.

Slow kettle cooked

Ingredients:

Beef, Water, Diced Tomatoes (Tomatoes, Tomato Juice, Salt, Citric Acid, And Calcium Chloride), Dark Red Kidney Beans (Dark Red Kidney Beans, Water, Sugar, Salt, Calcium Chloride, Disodium Edta [To Preserve Color]), Tomato Paste, Onions, Green Peppers, Celery, Salt, Flavorings And Coloring, Chili Powder Blend (Chili Pepper, Salt, Spices, Garlic), Brown Sugar, Garlic Powder.

Serving Suggestions:

- Top with shredded cheddar jack cheese, diced onions, jalapeños and sour cream.
- Serve over French fries, tortilla chips, or baked potato.
- Serve with freshly baked cornbread for a perfect meal.

Product Information:

Unit Weight: 4 x 4 pound

Count: 4 (8 x 8 oz. serving) or 256 oz. per case **Shelf Life:** 12 months frozen or 7 days refrigerated



Nutrition Facts Serving Size 1 Cup (245g) Servings Per Container About 30 Calories 280 Calories from Fat 140 % Daily Value* Total Fat 16g Saturated Fat 6q 31% Trans Fat 0g 17% Cholesterol 50mg Sodium 880mg 37% Total Carbohydrate 16g 5% Dietary Fiber 5g 18% Sugars 5g Protein 18g Vitamin A 15% Vitamin C 15% Calcium 4% Iron 15% *Percent Daily Values are based on a 2,000 calorie diet. Your daily values may be higher or lower depending on your calorie needs



Packaging Information:

Configuration: 4 x 4 pound bags per case

Case Net Weight: 16 lbs. Case Gross Weight: 17 lbs.

Case Dimensions: 13.12" L x 8.37" W x 5.62" D

Case Cube: 0.36

Palletization: 16 TI x 7 HI

Preparation Instructions:

Fully cooked, heat & serve.

For best results, defrost overnight under refrigeration.

To Heat: Place unopened bag in a large pot of simmering water and heat until 160°F is reached.

To Take Temperature: Carefully remove bag from water, lay flat. Place thermometer in the center and fold bag around thermometer.